

obi[®]



INSTRUCTIONS FOR USE

IFD-600-003-0, INSTRUCTIONS FOR USE, R5.3.1 (ECO#185)

APPLICABILITY

These Instructions for Use cover Obi Gen1, Obi Gen2, and Obi3 (Gen3). Certain features, utensils, controls, specifications, and operating instructions vary by generation. Generation-specific information is identified where applicable. Confirm the model shown on the device label before following generation-specific instructions.

Device version	UDI-DI/GTIN
Obi Gen1 and Gen2	00869174000106
Obi Gen3/Obi3	00869174000113

ABOUT DESIN LLC

DESIN LLC was founded in response to a personal family health tragedy in which DESIN LLC's founder witnessed his grandfather lose the ability to eat independently due to a degenerative neuromuscular disease. He painfully watched as his grandfather became depressed and especially apprehensive about his daily meals. After observing the same issue in nursing homes, and even with young patients, and unsatisfied with limited technology available on marketplace, DESIN has made it its mission since 2010 to lead the creation of a first of its kind medical device, Obi, the robot dining companion.

It is our sincere mission to provide an unparalleled experience and benefits to everyone in need of access to lasting independence.

DISCLAIMER

This Instructions for Use is a controlled document. Confirm that you are using the current revision available from DESIN LLC before relying on an electronically stored or reproduced copy. Changes are made periodically to the information herein; these changes will be incorporated into new revisions of this publication. DESIN LLC may make improvements and/or changes to the products and/or software programs described in this publication at any time.

If any information in this Instructions for Use is unclear or appears inconsistent, discontinue use until clarification is obtained from DESIN LLC or its authorized representative.

COPYRIGHT NOTICE

© 2016 DESIN LLC. All rights reserved.

This instructions for use is only to be copied or used according to terms of sale. Except as permitted by agreement, no part of this publication may be reproduced, stored in any retrieval system, or transmitted, in any form or by any means, electronic, mechanical, recording, or otherwise, without prior written consent of DESIN LLC.

This product uses software/firmware to which the company holds the rights.

Address comments to: info@meetobi.com

TRADEMARKS

® Obi is a registered trademark of DESIN LLC.

PATENTED

US 8,442,669 B2; US 8,761,922 B2; US 8,918,200 B2; US 9,101,238 B2; US 9,504,343 B2; US 9,008,826 B2; US 9,566,714 B2; US D740,071 S; US D696,325 S; US D740,075 S. Other patents pending.

NOTES

This page is intentionally blank.

NOTES

This page is intentionally blank.

TABLE OF CONTENTS

APPLICABILITY	2
ABOUT DESIN LLC.....	3
DISCLAIMER.....	3
COPYRIGHT NOTICE	3
TRADEMARKS.....	3
PATENTED.....	3
NOTES	4
TABLE OF CONTENTS.....	6
INTRODUCTION.....	8
IMPORTANT SAFEGUARDS.....	9
REGARDING THE LITHIUM-ION BATTERY	14
SYMBOLS	15
PACKAGE CONTENTS	18
INTENDED USE	19
INDICATIONS FOR USE / INTENDED PURPOSE.....	20
ABSOLUTE CONTRAINDICATIONS.....	20
RELATIVE CONTRAINDICATIONS/PRECAUTIONS.....	21
SIDE EFFECTS	21
INTENDED USE ENVIRONMENT	22
PRODUCT DESCRIPTION.....	23
OPERATING INSTRUCTIONS	24
FOOD PREPARATION.....	24
MICROWAVE.....	26
SET UP	27
UTENSIL SELECTION (OBI3 ONLY)	28
POWER ON	29
PATIENT POWER ON (OBI3 ONLY)	30
TEACH MODE.....	30
REWARD MODE (OBI3 ONLY)	32
DELIVERY COMMAND.....	32
CHOOSE COMMAND	33
RETURN COMMAND	33
FOOD REMOVAL & REPLACEMENT COMMAND	33
REPOSITION SEQUENCE	33
FOOD MIXING/DIPPING SEQUENCE.....	33

ONE SWITCH MODE	33
SLEEP MODE.....	33
COLLISION DETECTION	34
POWER OFF	34
BATTERY CHARGE	35
CLEAN & DISINFECT	36
<i>FOOD CONTACT PARTS.....</i>	<i>36</i>
<i>BASE, ARM & SWITCHES.....</i>	<i>37</i>
MAINTENANCE	38
STORAGE.....	38
TRANSPORT.....	39
SHIPPING.....	40
END OF LIFE.....	41
CONTACT INFORMATION.....	42
APPENDIX I: LIGHTING PATTERNS	43
APPENDIX II: TROUBLESHOOTING.....	44
APPENDIX III: TECHNICAL SPECIFICATIONS	45
APPENDIX IV: FCC COMPLIANCE STATEMENT	46
APPENDIX V: EU DECLARATION OF CONFORMITY AVAILABILITY	47
APPENDIX VI: EU FOOD CONTACT DECLARATION OF COMPLIANCE.....	48
APPENDIX VII: LIMITED WARRANTY.....	49
APPENDIX VIII: TERMS AND CONDITIONS	50

INTRODUCTION

Obi is a reusable robotic utensil intended to restore functional self-feeding—an essential activity of daily living.

Obi provides an eligible user the function to pick up a utensil, use suitable utensils, scoop food onto a utensil, bring food to the mouth, and deliver and drink fluids. The user operates Obi through its integrated controls or compatible accessibility switches. The delivery position can be established using Teach Mode and stored by the device to support consistent use and reduce the need for repeated setup.

Obi is intended to be used only by individuals who meet the criteria described in the “Indications for Use / Intended Purpose” section. Before Obi is prescribed or provided for home use, the prospective user should undergo the clinical assessment described in the “Intended Use” section. Refer to the “Contraindications,” “Precautions,” and “Important Safeguards” sections for additional information concerning appropriate user selection and safe operation.

About These Instructions for Use

These Instructions for Use contain the information necessary to set up, operate, clean, maintain, store, transport, and dispose of Obi safely and properly. They also describe the device’s intended use, indications, contraindications, precautions, technical specifications, and applicable safety warnings.

Read these Instructions for Use carefully before setting up or operating Obi. Understanding and following the instructions, warnings, and precautions will help promote safe and effective use. Keep this document in a convenient location for future reference.

An electronic version of these Instructions for Use is available at www.meetobi.com.

Regulatory, Safety, and Quality Conformity

Obi is designed, developed, manufactured, distributed and serviced under an ISO 13485 quality management system certified by EAGLE Registrations. Obi is also covered by an active SGS North America, Inc. safety certification listing (No. 801023), based on evaluation against applicable IEC 60601-1 and associated collateral-standard requirements, and bears the applicable SGS certification mark. It has been assessed for conformity with Regulation (EU) 2017/745 and bears the CE marking accordingly. The current EU Declaration of Conformity is available from DESIN LLC, its European authorized representative, or www.meetobi.com.

IMPORTANT SAFEGUARDS

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

When using electrical products, especially when children are present, basic safety precautions should always be followed, including the following:

READ ALL INSTRUCTIONS BEFORE USING

DANGER:

To reduce the risk of electrocution, burns and injury:

- DO NOT use this device in a manner inconsistent with its' intended use and operating instructions.
- DO NOT use Obi as a method of administering medication.
- Never use Obi while sleeping or drowsy.
- DO NOT use Obi while bathing.
- DO NOT place or store Obi where it can fall or be pulled into a tub or sink.
- DO NOT place in or drop Obi into water, other liquid or snow. Obi is water resistant only (IP21)
- DO NOT reach for Obi if it has fallen into water. Unplug immediately.
- Check the power supply cable for damage before each use.

WARNING:

To reduce the risk of electrocution, burns, fire, or injury to persons / device:

Power supply cable:

- Never operate Obi if: it has damaged cords or plugs; if it is not working properly; if it has been damaged or dropped into water. Return the product to a service center for examination and repair. Please see the "Contact Information" section.
- Keep the cords away from heated surfaces.
- DO NOT use the power supply cable in steamy, damp, or dusty areas.

- DO NOT modify the power supply cable in any way.
- DO NOT cover the power supply cable with clothing, blankets, or other materials.
- DO NOT touch the power supply cable with wet hands.
- DO NOT plug this device into a voltage system that is different from the voltage system specified on Obi or the Obi power supply cable.
- DO NOT place the power supply cable within any bodily orifice.
- DO NOT use or leave the power supply cable in direct sunlight.
- DO NOT touch the power supply cable for a long period of time.
- Only use the power supply cable provided to charge Obi:

Power Partners, PEAMD-72SF-13-1-BT1

- If the Obi power supply cable does not fit the wall outlet in your region, please contact the distributor for a different power supply or a recommended power adapter. Only use a power adapter recommended by the distributor or DESIGN LLC.
- Unplug the power supply cable during lightning storms or when unused for long periods of time.

Children:

- Close supervision is highly recommended when initially observing or evaluating a child's ability to operate Obi safely and successfully. Young children, because of their small fingers, are at greater risk of being pinched by the Arm.
- To avoid danger of suffocation, keep plastic bags away from babies and children.
- DO NOT use bags in cribs, beds, carriages, or playpens. The bags are not toys.
- CHOKING HAZARD – small parts. Not for children under 3 years. Use with adult supervision only.

Electromagnetic Energy:

- Be aware that adding accessories or components close to Obi's base, may make it more susceptible to electromagnetic interference (EMI) (There is no easy way to evaluate their effect on the overall immunity of the power wheelchair and/or your Obi).

- Be aware of nearby transmitters, such as radio or TV stations, and try to avoid coming close to them.
- In the presence of EMI, the following degradations may occur, and are free from unacceptable safety risk to the patient or user:
 - Motion disruption by arm falling with all visual indicators flashing, then powering off.
 - Movement from above one bowl to above another (choose sequence), without a user command.
 - Scoping movement into a bowl, followed by movement to the point of delivery (food delivery sequence), without a user command.
- Report all incidents of unintended shut down to the manufacturer or a local distributor and note whether there is a source of EMI nearby.

**Environment of Use:**

- DO NOT use where aerosol (spray) products are being used.
- DO NOT use in a pressurized oxygen rich environment, such as a hyperbaric oxygen chamber.
- DO NOT operate or store Obi outdoors or in extreme environments. Obi is for indoor use only.
- DO NOT place Obi on a wet surface.
- DO NOT place, store or operate Obi on an inclined surface.
- DO NOT place containers of liquid or any other item within 12 inches of the device to prevent interference with the arm during motion.

**Food Preparation:**

- Unplug Obi before filling with foods or liquids.
- Only fill the Plate's bowls with food and beverages.
- DO NOT place any component of Obi in an oven.

**Microwave Heating:**

- DO NOT microwave any other portion of Obi other than the Plate.
- When heating food, DO NOT use a power setting greater than 600 W.
- Do not microwave above 70°C (158°F).
- If any food in contact with the Plate reaches a temperature of 150°C (300°F) or greater, damage to the Plate is possible (blistering, discoloration, etc.).
- Microwave at your own risk of damage to the Plate.
- Do not cook food in the Plate; reheat only.

**Miscellaneous:**

- DO NOT carry Obi by the Arm.
- DO NOT place any part of Obi, other than the Plate, Placemat, and Spoons, in the dishwasher.
- DO NOT use solvents on Obi (example: turpentine, mineral spirits, gasoline, petroleum distillates, etc.). See Operating Instructions for proper cleaning.
- DO NOT cross or tangle wires connected to Obi to prevent trip hazards.
- Wash the Obi Spoon, Plate, and Placemat before their first use and after each subsequent use.
- To maximize Obi's lifespan, avoid continuously keeping the robot in motion for periods longer than one hour at a time.
- The product label with the serial number is located on the bottom (sitting surface) of Obi. See "Technical Specifications" for more information.

**Modifications:**

- Any modification or disassembly of the device may void the product's warranty.
- DO NOT modify Obi in any way as this may damage the device or lead to severe injury.
- Only use attachments recommended by the manufacturer.

- Modifications made by the user not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.
- DO NOT intentionally deform or misshape the Spoons. The product will no longer work.
- Never drop or insert any object into any opening.
- DO NOT push, pull, or force the Arm to move while Obi is powered on.
- DO NOT remove any labeling from Obi.
- DO NOT place anything on the Arm. This may damage the Arm.
- DO NOT crack or damage the housing of Obi. If the housing of Obi is cracked or damaged, DO NOT operate the device and contact the manufacturer immediately.
- Improper assembly of the Plate and Placemat on Obi may cause damage to the device and spillage of food.
- DO NOT drop Obi – this could degrade the product's performance. Handle with care and only as directed.
- Save the original packaging box! Obi is a delicate precision instrument, DESIGN LLC will not be responsible for the failure of custom packaging to adequately protect Obi from damage.

**Service:**

- DO NOT attempt to disassemble Obi, this will void the warranty.
- Only service the device with an authorized service representative.
- Return the product in the original packaging box.

SAVE THESE INSTRUCTIONS

REGARDING THE LITHIUM-ION BATTERY



For Customers in the U.S.A. and Europe:

THIS PRODUCT CONTAINS A 49Wh LITHIUM ION BATTERY.

DISPOSE OF LITHIUM-ION BATTERIES PROPERLY.

DO NOT attempt to disassemble Obi and remove or replace the battery pack.

Contact a local recycler, waste treatment facility, the distributor or DESIN LLC for further information on proper disposal. Please see the “Contact Information” section.



WARNING:

- The battery pack used in this device may present a fire or chemical burn hazard if mishandled.
- DO NOT dispose of batteries in a fire or trash incinerator.
- Use of a different battery may present a risk of fire or explosion.
- DO NOT disassemble, puncture, modify, drop, throw, or provide unnecessary stress to batteries.
- The battery pack used in this device can be recycled and a local waste treatment center may accept it.
- Dispose of batteries according to local laws and regulations in your region.
- Some regions prohibit the disposal of lithium-ion batteries in household and business trash.
- If you are not able to identify the applicable rules or waste treatment centers in your area, please contact a distributor or DESIN LLC.
- Special shipping requirements for lithium-ion batteries may apply.




CAUTION:

- Charge the battery at least once a year for infrequent use.
- Ensure the battery has at least a 40% nominal charge (2 solid LEDs) prior to storing for prolonged periods. This will help prevent a dead battery and increase battery life.

SYMBOLS

The following is a list of symbols that appear on the product, its packaging or labeling:

Symbol	Description
	ISO 15223-1, No. 5.4.4: Caution. Consult the instructions for use for important cautionary information such as warnings and precautions. Also used to highlight the fact that there are specific warnings or precautions associated with the medical device, which are not otherwise found on the label.
	ISO 15223-1, No. 5.1.1: Manufacturer.
	ISO 15223-1, No. 5.1.2: Authorized representative.
	ISO 15223-1, No. 5.1.3: Date of manufacture.
	ISO 15223-1, No. 5.1.10: Model number.
	ISO 15223-1, No. 5.3.1: Fragile, handle with care.
	ISO 15223-1, No. 5.3.4: Keep dry.
	ISO 15223-1, No. 5.3.7: Temperature limit. Indicates the temperature limits to which the medical device can be safely exposed.
	ISO 15223-1, No. 5.3.8: Humidity limitation. Indicates the range of humidity to which the medical device can be safely exposed.

	<i>ISO 15223-1, 5.4.3: Consult instructions for use.</i>
	<i>ISO 15223-1, No. 5.7.7: This device is a Medical Device (MD).</i>
	<i>IEC 60417-5031: Direct current.</i>
	<i>IEC 60417-5333: Type-BF Applied Part (device).</i>
	<i>ISO 7000, No. 0623: This way up.</i>
	<i>ISO 7010-M002: Refer to instruction manual – mandatory action.</i>
IP21	<i>Degree of ingress protection. Protected against access to hazardous parts with a finger, against solid foreign objects 12.5 mm in diameter and larger, and against vertically falling water drops</i>
	Indicates that Obi is covered by an active safety certification listing with <i>SGS North America, Inc.</i> and has been evaluated against the applicable medical electrical safety requirements in the IEC 60601-1 family and identified in the SGS certification record. The number 801023 identifies the applicable SGS listing.
	The letters " CE " abbreviate the French phrase <i>Conformité Européenne</i> , meaning, <i>European Conformity</i> . It signifies that the device complies with applicable directives and regulations of the European Union (see Declaration of Conformity).

	For users in the European Union: Compliance with the <i>Waste Electrical & Electronic Equipment (WEEE) Directive, 2012/19/EU, as amended.</i>
	This device complies with the European <i>Restriction of Hazardous Substances Directive (RoHS) Directive 2011/65/EU, as amended, including Commission Delegated Directive (EU) 2015/863.</i>
	This part is appropriate for food contact, with the exception of cooking food, and complies with regulation (EC) No 1935/2004 and (EU) No 10/2011, as amended.
	Microwave compatible: In the presence of microwave radiation, the Plate alone will not heat to an unsafe handling temperature, warp, melt, or damage the microwave oven.
	Dishwasher compatible: Within a dishwasher top rack environment, the part alone will not heat to an unsafe handling temperature, warp, melt, or damage the dishwasher.
	This part is safe to store in the refrigerator and freezer (Plate).
	Please recycle this part at the end of its life (Plate, Placemat, & Spoons).
	This part is BPA free (Plate, Placemat, & Spoons).

PACKAGE CONTENTS

The Obi package includes the following contents:

- Obi
- Plate(s)
- Placemat
- Utensils
- 2X Accessibility Switches
- Power supply cable & Power Supply
- Instructions for use
- Quick Start Guide



PLEASE SAVE OBI'S ORIGINAL FOAM AND PACKAGING MATERIALS!

Obi is a precisely calibrated device and custom packaging solutions may not adequately protect Obi from functional damage.

INTENDED USE

Obi is a reusable robotic utensil intended to restore functional self-feeding—an essential activity of daily living.



CAUTION:

- Before Obi is prescribed or provided for home use, a qualified healthcare professional knowledgeable about assistive feeding devices should evaluate the prospective user with Obi during an in-person clinical assessment. The assessment should determine whether the individual meets the indications for use, can use Obi safely and effectively in the intended home environment, and what training, assistance, supervision, positioning, or patient-specific precautions are required. Continued appropriateness should be reassessed whenever the user's physical, cognitive, behavioral, swallowing, positioning, caregiving, or environmental circumstances materially change.
- Obi does not eliminate any prescribed/pre-existing dietary restrictions.
- Although Obi can provide liquids, for some individuals, additional methods of hydration may be necessary to ensure adequate hydration.

INDICATIONS FOR USE / INTENDED PURPOSE

Obi is indicated for people 5 years and older (or younger, with approval of a qualified healthcare professional) living with severe congenital or acquired neurological, neuromuscular, or musculoskeletal impairment of the upper extremities affecting their ability to functionally self-feed. Obi users may have the following medical diseases/conditions:

- A. Amputees
- B. Amyotrophic Lateral Sclerosis (ALS)
- C. Arthrogryposis Multiplex Congenital (AMC)
- D. Cerebral Palsy
- E. Congenital Limb Differences
- F. Developmental Disabilities
- G. Essential Tremor
- H. Multiple Sclerosis
- I. Muscular Dystrophy
- J. Parkinson's Disease
- K. Spinal Cord Injury
- L. Other conditions impairing motor control of the upper extremities.

In addition, patients must have the cognitive ability to understand cause and effect (i.e. activate a switch to perform a function).

ABSOLUTE CONTRAINDICATIONS

Obi has no known absolute contraindications as there are no conditions under which the device must not be used because the risk of use clearly outweighs any possible benefit. However, there are several conditions that may carry risk and therefore caution should be considered via a benefit/risk assessment with a qualified healthcare professional (see "Precautions").

RELATIVE CONTRAINDICATIONS/PRECAUTIONS

A patient should NOT use Obi if any of the following conditions are present and of a severity that results in an unacceptable compromise to the benefit/risk ratio provided by the device. Where appropriate, special attention should be given to whether the following attributes pose an unacceptable risk to unsafe feeding or eating.

- A. **Cognitive/behavioral impairment:** Inability to understand or perform correct use of the device and/or at risk for combative/self-injurious behavior during meals.
- B. **Deviations from pre-existing dietary restriction:** No changes to food preparation should be made to use Obi that would directly contradict any current recommended food or liquid textures/consistencies prescribed by a healthcare professional.
- C. **Dysphagia:** Significant oral motor, breath, or swallowing impairment which could pose an elevated risk of choking/aspiration (unless found to be appropriate by a treating healthcare professional who has reviewed and discussed the benefit/risk assessment with the patient and /or caregiver). Any specific recommendations by the provider must be strictly followed.
- D. **Excessive movement:** Uncontrollable/forceful movements that prevent correct and/or safe use of the device.
- E. **Insufficient head control:** Insufficient head control to reliably position the Obi utensil near the mouth. Or, insufficient head and oral motor control to reliably remove content from the utensil.
- F. **Insufficient movement:** Unable to operate any compatible accessibility switches: Pressure Switch, Toggle Switch, Pillow Switch, Sip and Puff, etc.
- G. **Reclined positions:** Patients that cannot be positioned sufficiently upright for food/liquid consumption or maintain a sufficiently upright position throughout a meal (unless found to be appropriate by a treating healthcare professional).
- H. **Open mouth sores:** Any significant cuts, ulcerations, or oral bleeding.
- I. **Open skin sores:** Any open wound on an area of the body used to activate Accessibility Switch(es) or to contact Obi.

SIDE EFFECTS

No undesirable side effects attributable to the device have been identified. Known and foreseeable risks associated with eating, swallowing, food preparation, positioning, and device operation are addressed in the warnings and precautions section of this Instructions for Use.

INTENDED USE ENVIRONMENT

Obi is intended to be used in a variety of environments and locations including (but not limited to):

1. Residential and Commercial Facilities:

- a. Assisted living centers
- b. Homes
- c. Hospices
- d. Hospitals
- e. Long Term Care/Nursing homes
- f. Rehabilitation centers
- g. Restaurants
- h. Schools / Universities

2. Locations:

- a. Cafeterias
- b. Classrooms
- c. Dining Rooms
- d. Kitchens
- e. Living room
- f. Restaurants

3. Mounting: Use the device on a surface meeting all the following characteristics:

- a. Immobile, flat and dry.
- b. Cleaned or disinfected sufficiently for an eating environment and the user's needs.

4. Mobility: If the user is in a wheelchair at the time of use, ensure the wheelchair is positioned in the following manner:

- a. Fixed/locked in an appropriate location for use.
- b. Oriented in a sufficiently upright position promoting safe and effective swallowing behavior (see Precautions).

5. Hygienic:

- a. Clean and disinfect the Plate, Placemat and Spoons in accordance with the CLEANING & DISINFECTION sections of these instructions for use.
- b. Follow basic hygiene practices, such as vigorous hand washing with soap and water, before and after handling food and food contact parts.

6. Frequency of Use:


- a. Obi is intended for intermittent use; to be operated continuously for an hour interval or less.
- b. Between each meal, at least a one-hour power-off break is good practice. Sufficient power off breaks help maximize the life of Obi.

PRODUCT DESCRIPTION

Obi is a machine for individuals who cannot feed themselves. Obi's functionality mimics a human arm during eating. At least one off-the-shelf Accessibility Switch is required to eat with Obi. There are various types and styles of Accessibility Switches for purchase through DESIN LLC and other vendors. Obi easily accommodates these, allowing the user to:

- Select what food to eat (Choose Switch).
- Command when to eat (Deliver Switch).

The following image illustrates Obi's main parts and features:



OPERATING INSTRUCTIONS


Prior to first use, please charge the device for approximately half an hour.

Please refer to the **Quick Start Guide** for general operating instructions and added visuals. If needed, download an electronic copy of the **Quick Start Guide** at www.meetobi.com

FOOD PREPARATION

If prepared properly, nearly any food works with Obi (solids, purees or liquids). Obi is also compatible with certain dysphagia diets + *IDDSI*¹ diets of Regular, Easy to Chew (cut up to appropriate size) and below.


Solids: For foods in solid form, please ensure to cut or prepare the food into pieces smaller than 19 mm x 19 mm ($\frac{3}{4}$ " by $\frac{3}{4}$ " or about the size of a raspberry. Optimally, 6 mm x 6 mm ($\frac{1}{4}$ " by $\frac{1}{4}$ " or about the size of a pea works best.



Semi-Solids/Purees: Foods that are prone to become sticky or coagulate (oatmeal, cheesy foods, mashed potatoes, etc.) work best when moist/thinned as this promotes ease of separation into appropriate bite sizes.

¹<https://www.iddsi.org/standards/framework>

Food Insertion: After ensuring food is of an appropriate size or consistency, place the food in each bowl as shown. For the best results, fill each bowl to the Optimal Fill Line molded into the Plate. For illustrative purposes, the Optimal Fill Line position is indicated below in red.



Note1: The Plate (with or without food) may be safely stored in the refrigerator or freezer. With food, it's recommended not to exceed 24 hours.



CAUTION:

- **Significantly overfilling or underfilling any bowl from the Optimal Fill Line is likely to degrade Obi's ability to deliver food.**
- DO NOT cook food in the Plate. The Plate is for serving and warming food only.
- Always vigorously wash hands with soap and water prior to preparing food or handling food contact parts.
- DO NOT place excessively hot food or liquids in the Plate which could potentially burn the user / patient.
- Use care if placing the Plate under a heat lamp. Do not leave under the lamp for extended periods as the part could soften and warp and will not be covered under warranty.

MICROWAVE

Use the microwave to warm or reheat food or liquid to a desirable eating temperature.

The Plate is the only Obi component deemed “microwave safe”. When the Plate alone is in presence of microwave radiation, it will not heat to an unsafe handling temperature, warp, melt or damage the microwave oven.

When heating food in the Plate, use a low power setting (≤ 600 W).

Microwave in short intervals (≤ 1 minute), stirring frequently in-between each interval to avoid hot spots.

Note1: To maximize Plate life, warm or cook food in other glass or ceramic cookware.



CAUTION:


- DO NOT cook food in the Plate.
- DO NOT microwave any other portion of Obi other than the Plate.
- When heating food, DO NOT use a power setting greater than >600 W.
- DO NOT microwave past 70°C (158°F) – migration has not been evaluated above this limit.
- If any food in contact with the Plate reaches a temperature of 150°C (300°F) or greater, damage to the Plate is possible (blistering, discoloration, etc.).
- Microwave at your own risk of damage to the Plate. Microwave damage will not be covered under warranty.

SET UP

Place Obi on an immobile, flat and dry surface. Ensure there are no dinner glasses or other objects within approximately 25cm (12in) of any part of Obi's Base.

Next, to ensure the robotic arm can reach the user's mouth, confirm the patient is positioned within 40cm (16") from the front of Obi.


Ensure the Arm is positioned similarly to the image below; Move to the right of the four



bowls. Obi will not turn on unless it's approximately oriented as shown.

Next, Insert the desired accessibility switches into the Deliver Port and/or Choose Port of Obi until a click is felt and heard. Obi also operates with a single Accessibility Switch installed. See ONE SWITCH MODE section.

Next, insert the Placemat and Plate ensuring the keyed slots align. When inserted correctly, the plate sits in an immobile manner.



CAUTION:

- DO NOT place containers of liquid or any other item within 25cm (12in) of the device to prevent obstructions with the Arm when it's in motion.
- DO NOT place Obi on a wet surface.
- DO NOT place Obi on an inclined surface.

Note1: Obi's Switch Ports have 3.5mm barrel jack (headphone jack) connectors. Any wired, un-powered accessibility switch, less than 153cm (5ft) in length is recommended for use.

Note2: Velcro tape may be helpful to secure adaptive switches in place i.e., headrest, floor, tabletop, armrest etc.

Note3: It's possible to mount Obi to a wheelchair. Wheelchair specific table and tray attachments may be available to use with Obi. Please contact your wheelchair manufacturer to see if such an accessory is available.

UTENSIL SELECTION (Obi3 Only)

The following are basic guidelines to help select the most suitable utensil(s). Always defer to the advice of a qualified healthcare professional. Small Spoon, Spork and Large Spoon are not available for Obi Gen1 and Gen2.

SMALL	CLASSIC	SPORK	LARGE
			
Note: Images are not to scale.			

For children 12 years and younger: To promote safe use, it's recommended to start with the Small Spoon. Following experience with the device, increase the size as appropriate.

For others: To promote safe use, it's recommended to start with the Classic Spoon. Following experience with the device, increase or decrease the size as appropriate.

The Utensils are keyed and only attached in the correct direction.

**CAUTION:**


Please use caution when considering the suitability of the Spork and Large spoon. These utensils may not be suitable for children under 12, and for individuals living with significant spasticity or dysphagia.

POWER ON

Ensure the battery is sufficiently charged prior to use (at least 2 solid illuminated LEDs recommended).

Ensure the Arm is positioned similarly to the image below; to the right of the four bowls, then gently press the Power Button.

Obi will not turn on unless it's approximately oriented as shown below.



The button will illuminate a continuous white and Obi will wake to greet you. The Arm will then position itself above the first bowl (nearest the power button).

Note1: If the Teach Button flashes after pressing the Power Button, the Arm must be moved. Please ensure the Arm is in a similar area and orientation to the previous image.

Note 2: If the device indicates an Alternating Flash of Power and Teach Button, the device is too cold or too hot. Please allow the device to sit at room temperature for up to a half hour before powering on again.

**CAUTION:**

When the Placemat and Plate align incorrectly, they will remain mobile, and Obi will NOT function properly.

PATIENT POWER ON (Obi3 Only)

If deemed appropriate by a qualified healthcare professional, it is possible to power on Obi through the Choose or Delivery switch attachment. This feature is not available on Obi Gen1 or Gen2.


To use Patient Power On, first plug in the Obi power cable and ensure the Battery LEDs are illuminated. Next, press and hold either the Choose or Deliver Switch for approximately 2 seconds. Obi will start up and respond as if the Power Button has been pressed.

TEACH MODE

Teach Mode enables a caregiver to quickly tailor a custom food delivery location in front of the user's mouth. Before entering Teach Mode, ensure the patient is sitting upright, and that the table and chair are sufficiently fixed/locked, such that they will not roll away during a meal.

To enter Teach Mode, gently touch the Teach Button on the Base. When the button activates, it illuminates continuous white.

Grasp the Arm with a finger fully covering the Teach Icon. In response, the arm will release, and you will feel the weight of the arm – this is normal. While continuing to contact the Teach Icon, lift the Arm as shown in the images below.




When the Teach Icon is covered, the Arm will be free for the caregiver to move. The last two sections of the Arm will remain engaged and will keep the Spoon level when moving within the teach space. If resistance to movement is present, the thumb is probably not completely covering the Teach Icon.

Teach the spoon to deliver food in front of the user's mouth. It's recommended to first try positioning the spoon mid-chin level. This position results in favorable ergonomics for many patients. After achieving a desired location, release the Arm. Obi immediately holds position. If another change in position is necessary, simply grasp the Arm, with your thumb or finger covering the Teach Icon, and move as previously instructed. Repeat as many times as needed until a satisfactory position is achieved.

Once satisfied with the food delivery location, press the illuminated Teach Button on Obi's Base to store the position. Obi will remember the taught position until re-taught (irrespective of how many times powered on and off). If the Teach Button is flashing, the current position is invalid (see notes below) and cannot be stored.

Note1: Not all physical positions are valid for teaching a point of delivery. An approximate recommended teaching space (for the tip of the Spoon) is illustrated grey in the aerial view below. If the Arm has been moved outside the illustrated teaching space, the Teach Button on the Base will begin to flash, and Obi will shake / rotate his head side to side. This signal means: "please relocate the Arm".



Note2: If left unsupported in an extended invalid position, the Arm will fall, and the Teach Button will continue to flash. When this occurs, promptly re-teach the Arm as previously instructed. If left unsupported for several seconds, Obi will power off. If this happens, simply move the Arm to the right of the Plate as shown in the "Set Up" section and power on Obi again.

Note3: A continuously illuminated Teach Button with the head shaking / rotating side to side means, "the Arm is being handled incorrectly". A user is pulling or pushing the Arm without touching the Teach Icon.

Note4: During Teach Mode, the Arm should move easily. First, please ensure a finger is in contact with the Teach Icon.



WARNING:

- Choking Hazard: DO NOT Teach a delivery position such that the Spoon will penetrate the user's mouth.
- Choking Hazard: DO NOT Teach a delivery position to accommodate a reclined position unless it has been approved by a healthcare provider.



CAUTION:

- DO NOT force the Arm to move when teaching. If there are any issues, please see the "Contact Information" section of this instructions for use.
- DO NOT intentionally deform or misshape the Spoons as they will no longer scoop through the bowls correctly. Attempt to re-teach the arm instead.

REWARD MODE (Obi3 ONLY)

This feature provides positive reinforcement for favorable feeding behavior or self-feeding related accomplishments. When activated, Obi will flash the lights in a unique pattern and Obi will move through unique paths and motions for approximately 1 minute. Once the time has expired, Obi will automatically exit Reward Mode and return to the Ready Position above the previously used bowl. Stay at least 1 foot away from Obi's base when this mode is activated to avoid collision/interference with the arm or utensil.

To activate Reward Mode, press the Teach Button to enter Teach Mode. Without touching the arm, hold down either the Choose or Deliver Switch for approximately 2 seconds. Obi will immediately start moving and flashing the power, teach and battery LEDs (this is normal).

Note1: If the Arm Sensor is pressed after entering Teach Mode, the Reward Mode will not activate as it's assumed the intent of entering Teach Mode is to select a new food delivery location.

Note2: Once Reward Mode is activated, it will not stop until 1 minute has expired.

DELIVERY COMMAND

To capture and transport food and liquid to the taught location, activate the Deliver Switch:



Note1: After activating an Accessibility Switch, Obi will carry out the complete choose or deliver sequence. While the Arm is in motion, interruptions (besides collision) cannot

occur. That is, the Accessibility Switches and Teach Button will only operate when the Arm is stationary. However, the user may activate the Power Button at any time to Power Off the device.

CHOOSE COMMAND

To select a bowl, activate the Choose Switch:



RETURN COMMAND

After eating from the Spoon, either the Deliver Switch or the Choose Switch can be activated. Press the Deliver Switch to return the Arm to the same Bowl. Press the Choose Switch to return the Arm to a different Bowl.

FOOD REMOVAL & REPLACEMENT COMMAND

If the user initiates a Return Command immediately before or after the Spoon arrives at the food delivery location, the Arm will dispense the food back in the proper Bowl.

REPOSITION SEQUENCE

After several scoops from the same bowl, Obi will automatically use the Spoon to reposition food back to the center of the bowl.

FOOD MIXING/DIPPING SEQUENCE

This sequence enables mixing of food from one bowl to another, including dipping/dunking in sauces, etc. To conduct the food mixing sequence, simply command a Deliver Sequence of food to the point of delivery. Without consuming the food, command the “Choose” sequence and the arm will move the spoon to above the next bowl in a counter-clockwise direction. Keep commanding the choose sequence to move the spoon to the desired bowl. Command the Deliver sequence to mix/dip the food into the chosen bowl.

ONE SWITCH MODE

Obi can operate off one Accessibility Switch. The installed Accessibility Switch will act like a Deliver Switch. That is, Obi will automate the choosing of bowls (moving one bowl position counterclockwise).

SLEEP MODE

After two minutes of inactivity, the Arm will move to the Storage Position to save battery power. In this position, the Arm remains stationary, and the power button will stay illuminated. During this time, if the device senses any command, Obi will immediately

wake up and continue normal operation. If the half hour elapses without any command, Obi will automatically power off.

COLLISION DETECTION

If an object interferes with the moving Arm with an unexpected force (including food obstructions during a scoop), Obi will stop and slightly back away. During this situation, the Power Button will flash. If a collision occurs with an external object, remove the obstructing object. If a collision occurred while scooping, simply clear the collision by pressing either switch and Obi will continue the intended path.

Note1: If in-bowl collisions occur repeatedly, please check food is prepared in accordance with the Food Preparation section of this manual.

POWER OFF

To turn off Obi, gently touch the Power Button continuously for one second. Obi will proceed to the Storage Position and turn off.

Note1: If the power off command initiates while the Arm is in motion, the Power Button will repeatedly double flash. Obi will proceed to the Storage Position and turn off once present movement is finished.

BATTERY CHARGE

Prior to first time use, please charge the device for approximately half an hour.

If attempting to charge a used device with a fully depleted and unresponsive battery, please see the TROUBLESHOOTING section.

When only one Battery Indicator Light is illuminated and flashing, Obi needs to charge. The single light will begin flashing quickly when Obi only has a few minutes of power left.

If Obi does run out of power during operation, the Arm will automatically move to the Storage Position and Obi will power off.

To Charge in the US:

Connect the power supply cable to a 120V, 60hz socket and to the Obi Charge Port, located under the device Arm. See the “Product Description” section for a visual.

To Charge in the European Union:

Connect the power supply cable to a 230V, 50hz socket and to the Obi Charge Port, located under the device Arm. See the “Product Description” section for a visual.

When the illuminated Battery Indicator Light flashes slowly, the battery is successfully charging. As the battery charges from empty to full, the five Battery Indicator Lights will proportionately illuminate a continuous white. For example, if three of the five Battery Indicator Lights are illuminated a continuous white, the battery is sixty percent charged.

Once all five Battery Indicator Lights are illuminated, and no flashing occurs, Obi has a full charge. Disconnect the power supply cable from Obi and the wall socket.

NOTE 1: Obi does not need a full charge before disconnecting the power supply cable; partial charging is acceptable.



CAUTION:

- Only use the Power Supply provided with Obi: *PEAMD72-SF-13-1-BTI*
- If the Obi Power supply cable does not fit the wall outlet in your region, please contact the distributor for a different power cord adapter. Only use an AC power adapter recommended by DESIN LLC.

CLEAN & DISINFECT

The following instructions are recommendations for cleaning, disinfecting and sanitizing Obi's components. However, it is ultimately the responsibility of the customer/user to ensure the device is appropriately cleaned, disinfected or sanitized prior to use. For institutional customers/distributors facilitating reuse devices, this includes compliance with applicable laws related to food safety, food handling, personal hygiene and others.

FOOD CONTACT PARTS

Component	Home/Individual User	Institutional/Multiple Users
Plate, Placemat, & Utensils	A. Frequency: after each meal B. Method(s): - STEP 1: Rinse, and if necessary, scrub, components in the sink to remove bulk food particles. - STEP 2: Place parts on the top rack of a dishwasher (NSF/ANSI184 certified dishwasher recommended, or equivalent). Configure to sanitize setting (if applicable) - STEP 2 (Alternate): Scrub all surfaces with dish soap & warm water. Rinse.	A. Frequency: after each meal B. Method(s): - STEP 1: Rinse, and if necessary, scrub, components in the sink to remove bulk food particles. - STEP2: Sanitize in a commercial dishwasher meeting the NSF/ANSI3 standard (or equivalent).



CAUTION:

1. It is best practice for each new user of Obi to obtain their own set of Obi dishware.
2. If cleaning within an industrial dishwasher, ensure the temperature of the dishwasher does not exceed 100°C. Exceeding this temperature could degrade the parts.
3. When attempting a new method of cleaning or disinfection, first attempt in an inconspicuous area.
4. DO NOT use food contact parts with visible signs of degradation.
5. DO NOT rub any components of the device with highly abrasive materials such as metal scrub brushes.
6. DO NOT use disinfecting wipes to clean food contact parts.
7. DO NOT use food contact parts for which the level of cleanliness is unknown, such as: parts which have fallen on the floor, parts handled by people, or parts which have been left sitting out.


BASE, ARM & SWITCHES

Component	Home/Individual User	Institutional/Multiple Users
Base & Arm	A. Frequency: after each meal B. Method(s): - STEP 1: Saturate a clean paper towel or cotton cloth with warm water. Wring out thoroughly. - STEP 2: Wipe down the device to remove errant food particles. - STEP 3 (optional): Obtain the appropriate compatible cleaner and use in accordance with its labeling. See the "Cleaner & Disinfectant Compatibility" table.	A. Frequency: after each meal B. Method(s): - STEP 1: Saturate a clean paper towel or cotton cloth with warm water. Wring out thoroughly. - STEP 2: Wipe down the device to remove errant food particles. - STEP 3: Obtain the appropriate compatible disinfectant and use in accordance with its labeling. See the "Cleaner & Disinfectant Compatibility" table.
Accessibility Switches	Clean per the manufacturer's specifications.	Disinfect per the manufacturer's specifications.



CAUTION:

- Common household disinfecting wipes, such as *Clorox Disinfecting Wipes*, or *Lysol Disinfecting Wipes* are NOT COMPATIBLE with Obi as they may contain Ammonium Chloride (see "Cleaner & Disinfectant Compatibility").
- When attempting a new method of cleaning or disinfection, first attempt in an inconspicuous area.
- The product has not been validated with sterilization techniques such as steaming. Sterilize at your own risk of damage.
- DO NOT spray the Base or Arm directly with any liquid.
- DO NOT place Obi in the dishwasher, or under a faucet of running water.
- DO NOT rub any components of the device with highly abrasive materials as this may damage the plastics.
- DO NOT use solvents such as turpentine, mineral spirits, gasoline, or petroleum distillates.



CAUTION: DO NOT USE
Ammonium Chloride based
cleaners and others (see below)

Cleaner & Disinfectant Compatibility

Cleaning/disinfectant products containing the following chemicals marked "Compatible" should be used when cleaning the Obi Arm or Base. Do NOT use products containing chemicals in the "Avoid" column as this may damage the enclosure material (cracking, crazing, discoloration, etc). This is not a comprehensive list of compatible cleaners and disinfectants.

Cleaning/Disinfectant	Compatible	Avoid	Recommended Brands (US only)
Vinegar & Water (1:1 ratio)	✓		Standard item. No specific brand.
Dish Soap & Water (1:10 ratio)	✓		Standard item. No specific brand.
Isopropyl Alcohol (70.0%/90.0%)	✓		Standard item. No specific brand.
Sodium Hypochlorite (Bleach) (0.5%/1:10)	✓		<i>Clorox Healthcare, Bleach</i>
Hydrogen Peroxide, (1.4-3.0%)	✓		<i>Clorox Healthcare, Hydrogen Peroxide</i>
Acetone		X	NA
Ammonia		X	NA
Ammonium Chloride		X	NA
Formaldehyde		X	NA
Potassium Hydroxide		X	NA
Sodium Hydroxide		X	NA

MAINTENANCE

Part(s)	Activity / Method	Conditions / Frequency
Utensils (Gen3 only)	Replace when A or B have occurred. For heavy use, annual replacement recommended.	A. Approximately 1,000 dishwasher cycles. B. The parts have visual signs of degradation, such as: warp, cracks, crazing or discoloration.
Plate, Placemat and Spoons	Replace when A or B have occurred. For heavy use, annual replacement recommended.	A. Approximately 400 dishwasher cycles B. The parts have visual signs of degradation, such as, but not limited to: warp, blisters, fissures, cracks or discoloration.
Battery	Charge to 40-60% (2-3 SOLID LEDs)	At least once every 6 months, or twice per year (2X/Yr) charge for ~ 20 minutes to avoid a dead or degraded battery.
Barrel Jack Connectors & Accessibility Switch Connections	- STEP 1: Direct compressed, dry air into each barrel jack port to remove any bulk trapped debris. - STEP 2: Use a soft cloth to gently wipe clean the barrel jack pins of the Accessibility Switches.	At least once per year (1/Yr) recommended. If or when the Accessibility Switches become intermittently unresponsive; as needed.



CAUTION:

- DO NOT use food contact parts with visible signs of degradation.
- There is no "home serviceable" part inside the product. DO NOT open.
- If any fissure, crack or crevice traps food and removal is not achievable with the CLEANING AND DISINFECTION procedures of this instruction for use, please contact an authorized service representative.

STORAGE

Environmental Conditions for Storage:	-20° to 60°C, 15 to 93% humidity (short duration)
	12° to 30°C, 40 to 55% humidity (recommended)

Prior to long duration storage, ensure the battery has between 40-60% charge.

Store the Plate, Placemat and Spoons in a location where there is a high confidence in an appropriate level of cleanliness for clean dishware.

Note1: Keep the Obi package, as it makes a great storage container.




CAUTION:

- DO NOT store Obi near stovetops, fireplaces, and radiant heaters.
- DO NOT store Obi near a floor drain, or excessive source of humidity or moisture.

TRANSPORT

For ease of carrying, place the Arm in the Transport Position as shown below (with or without the Spoon, Plate and Placemat).



Because of Obi's curved contours, carrying is convenient in a multitude of positions other than the one pictured above.

Obi MEDICAL EQUIPMENT BAG (BACKPACK): The Obi medical equipment bag is a purpose built, reusable, cleanable (wipeable), backpack promoting effective portability, protection from environmental damage, protection from environmental contamination, and organization of Obi and its' components (switches, charger, plate, utensils, etc). The backpack is designed specifically to hold Obi (all generations), its' accessories, and a water bottle/thermos. Additionally, its size is within the 53cmX36cmX23cm (21"X14"X9") carry-on size limits posted by most airlines and contains straps and clips for effective wheelchair attachment. The Obi backpack must be used with the same foam that comes standard with Obi.

The Obi backpack must be carried by an adult or attached to a wheelchair; not suitable to be carried by children. Always carry the backpack with both straps to help avoid back strain.

Inspect, and if necessary, clean after each use (refer to wash label on backpack).



CAUTION: To avoid damage to the Obi MEB/device and potential injury to the user:

- DO NOT use the Obi backpack as checked luggage when flying.
- DO NOT have a child carry the Obi backpack
- DO NOT insert Obi without using the proper Obi foam.
- DO NOT use any other brand of bag/backpack for transporting Obi as it may result in damage.
- Use caution if a pre-existing back/spinal condition exists – a loaded bag may aggravate the issue.

- Tripping hazard, store away from areas with foot traffic.
- DO NOT store outside in harsh environments for long periods of time.
- Avoid dropping the backpack with Obi inside.

SHIPPING

When shipping Obi by ground, sea, or air please abide by the following:

Environmental Conditions for Transport:	-20° to 60°C, up to 93% humidity (short duration)
	12° to 30°C, 40 to 55% humidity (recommended)

Obi should ship in its original packaging box/backpack with an over-box



CAUTION:

- Obi is a precisely calibrated device and custom packaging solutions may not adequately protect Obi from functional damage.
- Please notify the shipper that Obi is a fragile, delicate instrument that necessitates handling with care.
- Obi encloses a 49 Wh lithium-ion battery. Special shipping requirements for lithium-ion batteries may apply.
- Always include a desiccant bag (for at least 1.50 cubic feet) for protection against humidity and condensation.

END OF LIFE

Private households - Information on Disposal for Users of WEEE



This symbol on the product and / or accompanying documents means that used electrical and electronic equipment (WEEE) should not mix with general household waste. For proper treatment, recovery, and recycling, please take this product to designated collection points where it will be accepted free of charge.

Alternatively, in some countries, you may be able to return your products to your local retailer upon purchase of an equivalent new product.

Disposing of this product correctly will help save valuable resources and prevent any potential negative effects on human health and the environment, which could otherwise arise from inappropriate waste handling.

Please contact your local authority for further details of your nearest designated collection point.

Penalties may be applicable for incorrect disposal of this waste, in accordance with your national legislation.

For professional users in the European Union:

If you wish to discard electrical and electronic equipment (EEE), please contact your dealer or supplier for further information. Please see the "Contact Information" section.

For disposal in countries outside of the European Union:

This symbol is only valid in the European Union (EU). If you wish to discard this product, please contact your local authorities or dealer, and ask for the correct method of disposal. Please see the "Contact Information" section.

EU Waste Electrical and Electronic Equipment (WEEE) Directive:

In August of 2005, the European Union (EU) implemented the EU WEEE Directive 2002/96/EC and later the WEEE Recast Directive 2012/19/EU requiring Producers of electronic and electrical equipment (EEE) to manage and finance the collection, reuse, recycling and to appropriately treat WEEE that the Producer places on the EU market after August 13, 2005. The goal of this directive is to minimize the volume of electrical and electronic waste disposal and to encourage re-use and recycling at the end of life.

DESIN LLC has met its national obligations to the EU WEEE Directive by registering in those countries to which DESIN LLC is an importer. DESIN LLC has also elected to join WEEE Compliance Schemes in some countries to help manage customer returns at end-of-life.


If you have purchased DESIN LLC -branded electrical or electronic products in the EU and are intending to discard these products at the end of their useful life, please DO NOT dispose of them with your other household or municipal waste. DESIN LLC has labeled its branded electronic products with the WEEE Symbol to alert our customers that products bearing this label should not be disposed of in a landfill or with municipal or household waste in the EU. Instead, please be aware that DESIN LLC is making a return and collection system available to you, free of transportation and recycling costs, for discarding these products.

CONTACT INFORMATION

If any serious safety incident has occurred in relation to Obi, report the incident to the manufacturer and the Food and Drug Administration (FDA). For users in Europe, report to the competent authority of the Member State in which the user and /or patient is established.

For Users in the United States:

Designed, manufactured, and distributed by:



DESIGN LLC (pronounced "Design")

SRN: US-MF-000008647

7018 AC Skinner Parkway, Suite 270

Jacksonville, FL 32256

Tel: 1-844-HelpObi (435-7624)

E-mail: service@meetobi.com

Internet: www.meetobi.com

In Europe:

Focal Meditech BV

SRN: NL-AR-000000549

Address: Droogdokkeneiland 19

5026SP Tilburg

Netherlands

Tel: +31 (0) 13-533 3103


Email: info@focalmeditech.nl

Website: www.focalmeditech.nl



This list does not include the contact information for all importers and distributors of Obi.

APPENDIX I: LIGHTING PATTERNS

NORMAL OPERATION		
Lighting Pattern		Meaning
	Solid Power Button	Obi is powered on and running normally.
	Flashing (slow) Battery Point LED	Obi is actively charging the battery. Obi could be on or off.
	Flashing (fast) Bottom Battery Point LED	Obi has a critically low battery. The device will soon automatically return to the sleep position and power-off.
	Double Flashing Power Button	Obi is in the process of powering down.
	Solid Power and Teach Button	Obi is in Teach Mode; Delivery and Choose commands will not respond.
	Solid Power Button, Flashing Teach Button	The arm is out of position. If attempting to power on the device, please move the arm to a valid start up position, see the POWER ON section. If in Teach Mode, please move the arm to a different position using the Teach Icon.
FAULT/ERROR CONDITION		
If these errors occur frequently, please contact an authorized service representative.		
Lighting Pattern		Meaning
	Flashing, All Battery LEDs Simultaneously	Shutdown Error, Servo Voltage: An electrical anomaly may have been detected by the battery, Power Supply or Servo.
	Flashing, Downward Moving Battery LED	Shutdown Error, Servo Over-Load. Too much force may have been applied to the arm while at a stationary position (not moving).
	Flashing, Power Button	Collision: Normal movement has been stopped due to an overload sensed by one or more servos. An overload is likely due to interference with an environmental object.
	Flashing, Alternating Power and Teach Button	Shutdown Error, Extreme Temperature: Obi is outside the acceptable operational temperature window. Obi was probably recently removed from an extreme environment such as a car in winter/summer.
	Flashing, All LEDs Simultaneously	Shutdown Error, Servo Communication. A signal anomaly may have been detected on the communication line.

APPENDIX II: TROUBLESHOOTING

If the applicable issue is not listed below, or if the problem persists after troubleshooting, please contact an Obi authorized service representative with the serial number and date of purchase. Some common problems and their solutions are as follows:

SYMPTOM	REMEDY
Battery – Does Not Charge	a. Ensure the correct power supply is being used (PEAMD-72SF-13-1-BT1) b. Confirm that the power supply and power supply cable are firmly connected together. c. If Obi has been unused for a long duration, the battery may have entered a protective standby mode. Please attempt the following to initiate a recovery charge: i. Insert the live power supply cable into the Obi charge port. ii. Remove the power supply cable from Obi then quickly re-connect it (within 3s). The bottom LED should start flashing. iii. Allow approximately 2 hours to charge. In this event, it is normal for Obi to take longer than normal to recharge.
Arm – Does Not Startup	a. Check to ensure the Arm is in a valid start up position (please see POWER ON section). If the Arm is in an invalid location, the Teach Button will flash, and movement will not start. b. Attempt a few slightly different arm locations to the right of the bowls.
Scooping / Food Delivery Issue	To ensure an optimal experience, food size, quantity and consistency should be reasonably workable by a spoon (see FOOD PREPARATION section): a. Size: Ensure food is prepared within the correct size: Optimally, ¼” by ¼” (6mm X 6mm), or about the size of a pea, works best. However, up to ¾” by ¾” (19mm X 19mm), or about the size of a raspberry, should also work. b. Quantity: Ensure food is not filled above the max fill line in each bowl. c. Consistency: Ensure food is appropriately un-sticky and can be easily separated by a spoon scoop.
Choose/Delivery Command - Unresponsive	a. Ensure both accessibility switches are fully and firmly inserted. b. Relocate the problem accessibility switch to the alternate switch port to determine whether the connector or the accessibility switch is the problem. c. Consult the MAINTENANCE section of this document and perform the listed activities for the barrel jack connectors.
Collisions	As appropriate, ensure: a. the spoon is not physically deformed from its proper shape. b. food bowls are free of any large/hard obstruction. c. there are no surrounding objects within Obi’s space of movement. d. Obi comes to a complete stop at the delivery point before taking a bite. If the problem persists – Obi is likely in need of calibration. Please contact an authorized service representative.

APPENDIX III: TECHNICAL SPECIFICATIONS

Accessibility Switch Connector	3.5mm male, mono barrel jack (0.138") 
Adhesive - Arm	Epoxy (Medical Grade; ISO-10993)
Applied Part (entire device)	Type-BF
Base Mounting Surface Material	Thermoplastic elastomer (Latex & nat. rubber free)
Battery Life per Charge	2-5 hours
Battery Type	Lithium ion (49Wh / 3.2Ah)
Degree of Ingress Protection	IP21—protected against access to hazardous parts with a finger, against solid foreign objects 12.5 mm in diameter and larger, and against vertically falling water drops
Dimensions (Transport Position)	30.5cm X 43.2cm X 15.2cm (12" X 17" X 6")
Dishwasher Safe Components	Plate, Placemat, Large Spoon and Small Spoon.
Enclosure Material	PC/ABS, UL94V-0 flame class rating
Food Bowl Volume Capacity	118mL (4 fl. Oz.) / Bowl (up to max fill line)
Food Safe Components	Plate, Placemat, Large Spoon and Small Spoon.
Food Safe Materials	Polypropylene & Copolyester
Lowest Food Delivery Height	5.0cm (2.0in) (below the mounting surface)
Maximum Food Delivery Height	38.1cm (15.3in) (above the mounting surface)
Maximum Horizontal Reach	420mm (16.5") (from the front of Obi)
Maximum Motor Speed	20 rpm
Medical Device Classification	European Union: Class I, active, non-sterile, non-measuring medical device under Regulation (EU) 2017/745. United States: Class I, 510(k) exempt, daily activity assist device under 21 CFR 890.5050; FDA product code ILC.
Microwave Safe Components	Plate
Environmental Operating Conditions	12° to 40°C, 15% to 93% humidity (short duration)
	12° to 30°C, 40% to 55% humidity (recommended)
	700 hPa to 1060 hPa atmospheric pressure
Peak Power	68 W
Recyclable components	Plate (#5), Placemat (#5), Large Spoon (#5), Small Spoon (#5), Lithium Ion Battery (RRC2054). Other components may also qualify. Talk to your local waste treatment center.
Safety Certification	SGS North America, Inc., Listing No. 801023.
Safety Evaluation, Electrical	Evaluated against the applicable requirements of IEC 60601-1 and associated collateral and normative standards.

Serial Number	The number is located on the bottom surface of the Base, as shown (exact appearance may differ by model): 
UDI-DI (GTIN issued by GS1)	Obi Gen1/Gen2: 00869174000106 Obi3 (Gen3): 00869174000113
Weight – Obi (no food)	3.49kg (7.7lbs)

APPENDIX IV: FCC COMPLIANCE STATEMENT

FOR UNITED STATES USERS:

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

CAUTION:

- Modifications made by the user not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.
- Operating Obi while the power supply is plugged in and actively charging may prevent compliance with 47CFR15 Class B for radiated emissions.

APPENDIX V: EU DECLARATION OF CONFORMITY AVAILABILITY

The current EU Declaration of Conformity applicable to this device is available upon request from DESIN LLC or its European authorized representative. It is also available at: www.meetobi.com/.

Obi has been assessed for conformity with Regulation (EU) 2017/745 on medical devices and bears the CE marking accordingly.

APPENDIX VI: EU FOOD CONTACT DECLARATION OF COMPLIANCE

EU Food Contact Declaration of Compliance				
Manufacturer		DESIN LLC 7018 AC Skinner Parkway, Suite 270 Jacksonville, FL 32256		
Product covered by this declaration:		IFD-500-000, Obi (all generations)		
Parts covered by this declaration		Obi3, Spork (IFD-400-108) Obi3, Small Spoon (IFD-400-102) Obi3, Large Spoon (IFD-400-101) Obi3, Classic Spoon (IFD-400-100) Utensil LG (IFD-400-062) Plate (IFD-400-027) Utensil SM (IFD-400-025)		
Date of declaration (mm/dd/yyyy):		08/03/2026 – addition of Obi3 utensils (supersedes 04/20/2020 initial declaration)		
Declaration of compliance with:				
• (EC) No. 1935/2004 (as amended) • (EC) No. 10/2011 (as amended) • (EC) No. 2023/2006 (as amended)				
The monomer(s), additive(s) and processing aids used in the manufacture of the listed polypropylene parts are authorized by the Regulation.				
Information about the compliance of substances used that are subject to any restriction or specification				
Compliance with Overall Migration Limits (OML)		Overall migration is less than 1mg/dm ² (limit is 10mg/dm ²) under standard testing conditions laid down in Regulation (EC) No. 10/2011 (see Note#1).		
Parts	Resin Ingredient	Substance	Specific Migration Limits (SMLs)	Test results or calculated estimate (see Note#1)
IFD-400-108 IFD-400-102 IFD-400-101 IFD-400-100	Monomer	TMCD	5mg/kg	<5mg/kg
IFD-400-062 IFD-400-025	Additive	Aluminum	1mg/kg	<1mg/kg
	Additive	Proprietary1	6mg/kg	
	Additive	Proprietary2	1.2mg/kg	
IFD-400-027	Polypropylene Homopolymer	DSTDP, CAS#693-36-7	5mg/kg	
	DSTDP Stabilizer	DSTDP, CAS#693-36-7	5mg/kg	
Information about the use of “dual use” additives in the material				
IFD-400-027	E470a, Calcium salts of fatty acids			
Conditions of use:				
Food Types		All		
Time & Temperature Limits		While in contact with foods:	70°C for 2 hours (do not exceed)	
		While stored in contact with food:	24 hours, <4°C (recommended)	
Information used to establish the compliance of the material or article				
The ratio of food contact surface area to volume	IFD-400-062	107.83(cm ²): 11.52(cm ³)		
	IFD-400-027	1599.02(cm ²): 188.51(cm ³)		
Functional Barrier				
There is no functional barrier present.				
Note#1: Further information is available with a non-disclosure agreement in place. Reference the following certificate of analysis: <i>DESIN LLC 20-75343 CoA ver.01</i> . Additionally, a comprehensive food contact report is available: <i>IFD-DES-024, Obi Food Contact Report</i> .				

APPENDIX VII: LIMITED WARRANTY

DESIN LLC ("DESIN") guarantees, subject to the terms below, that if your Obi is verified to be defective in material or workmanship at the time of delivery, DESIN will, at its discretion, either: (i) repair Obi using new or refurbished parts or (ii) replace Obi with a new or refurbished device.

Terms: United States -the limited warranty expires 12 months from the date of shipment or after the first 400 hours of usage, whichever occurs first. Europe - for warranty information, please contact your distributor.

Limitations: This Limited Warranty excludes Spoons, Plates, Placemats, and switches. Spoons, Plates, Placemats, and switches are non-returnable once used. Only manufacturing defects of materials or workmanship are covered when Obi is used properly and maintained in accordance to all included instructions. This warranty is void if Obi is improperly handled, dropped, put in a microwave, or exposed to excessive liquid(s), etc. In addition, any attempt(s) to open or dismantle Obi, including the Arm or Base assembly, will void this warranty. The warranty is limited to the initial purchaser of the Obi. The term of warranty identified above also operates to limit any implied warranties of merchantability or fitness under State law. Some States DO NOT allow limitations on how long an implied warranty lasts, so the above limitation may not apply to you. DESIN will not be liable for any incidental or consequential damages as a result of customer use. Some States DO NOT allow the exclusion or limitation of incidental or consequential damages, so the above limitation or exclusion may not apply to you.

Instructions: In order to obtain warranty service, deliver the Obi in either its original packaging or packaging with an equal degree of protection to the service center designated by DESIN or an authorized representative. An intact serial number or UPC is required for warranty claims. An issued RMA number may also be required by DESIN or an authorized distributor. The customer is responsible for all shipping and handling charges for the return of the Obi for warranty service, unless otherwise notified. For more information on how to obtain warranty service on your device in the US or Europe, please see "Contact Information", or contact your distributor.

Other Terms and Conditions/Legal Rights: This warranty is issued pursuant to the terms and conditions of sale entered into when you purchased your Obi, the provisions of which govern this Limited Warranty. This Limited Warranty gives you specific legal rights, and you may also have other rights which vary from State to State.

Return Policy: An unopened, unused or lightly used Obi generally may be returned within 14 days of the shipping date for a refund. A restocking fee may apply.

Please contact 1-844-HelpObi or Service@MeetObi.com for return details.

APPENDIX VIII: TERMS AND CONDITIONS

1. Acceptance: All Purchase Orders (as defined below) or other proposed agreements for the sale of goods ("Goods") by DESIN LLC ("DESIN") are subject to acceptance, in writing, by DESIN and are not binding on DESIN unless so accepted. These Standard Terms and Conditions ("Standard Terms of Sale") are applicable to all purchase orders, quotations, schedules or ancillary agreements concerning the sale of Goods, of any nature (collectively, the "Purchase Order"), and are the only terms and conditions that apply to the sale of DESIN's Goods, except those specifically set forth in any Purchase Order and accepted in writing by DESIN. Any term or condition or standard of performance different from or in addition to these Standard Terms of Sale, whether set forth on the Purchase Order or otherwise proposed by Customer, must be agreed to, in writing, by DESIN in advance, and in the event of any conflict or inconsistency, between these Standard Terms of Sale and any Purchase Order or other document proposed by Customer, these Standard Terms and Conditions will control. Customer agrees that Customer's submission of a Purchase Order to DESIN or Customer's acceptance of any Goods provided by DESIN shall constitute Customer's irrevocable acceptance of these Standard Terms of Sale.

2. Prices: Prices for the Goods set forth in the Purchase Order, unless otherwise stated in writing, DO NOT include taxes, as set forth in these Standard Terms and Conditions, transportation or shipping charges, special packaging or other services unless otherwise set forth in the Purchase Order. Prices are all quoted in denominations of United States' Dollars. The price for custom items not otherwise available for sale will be quoted by DESIN on a job-by-job basis. DESIN reserves the right to correct any error in prices.

3. Terms of Payment: Unless otherwise specified in writing, payment shall be received by DESIN prior to any shipment of any Goods. Payment shall be made in denominations of United States currency.

4. Customer Inspection: Upon Customer's receipt of shipment, Customer shall immediately inspect the Goods. Unless Customer provides DESIN with written notice of any claim for shortages of or patent defects in the Goods within five (5) days after receipt of shipment, such Goods shall be deemed finally inspected, checked and accepted by Customer and in full compliance with these Standard Terms of Sale.

5. Delivery/Packaging:

- a. DESIN shall use reasonable commercial efforts to make delivery of Goods or perform services within the time requested on Customer's Purchase Order. DESIN shall use its own discretion in choice of carrier and method of packing. DESIN has no obligation to obtain insurance for Customer covering Goods in transit to Customer.
- b. Prices stated are based on DESIN's standard packaging. Packaging will be standard commercial packaging and acceptable to the commercial carrier.

6. Cancellation and Returns: Only unopened and unused Goods can be returned to DESIN for a refund of the Customer's purchase price, provided Customer returns such Goods to DESIN within 14 days after the date of shipment to Customer. All returns must be postmarked by the fourteenth day. Customer is responsible for all shipping costs to return any Goods. Prior to issuing any refund, DESIN reserves the right to inspect all returned Goods to determine whether such Goods remain unopened and unused. No other returns of Goods shall be accepted without the prior written approval of DESIN. DESIN may exclude special orders from this 14-day return policy if done so in writing prior to accepting Customer's order.

7. Force Majeure: DESIN will not be liable or otherwise responsible for any damage, loss, fault, or expenses arising out of delays in manufacturing, shipment or other non-performance of any Purchase Order caused or imposed by: strikes, fires, disasters, riots, acts of nature; intervention of government, war or threat of war, acts of terrorism, conditions similar to war, sanctions, blockades, embargoes; acts of Customer; shortages of labor, fuel, power, materials, supplies, transportation, or manufacturing facilities; governmental action, subcontractor delay or any other cause, condition or circumstance beyond DESIN's reasonable control (collectively, "Force Majeure"). If there is a delay or nonperformance due to Force Majeure, then DESIN may, at its option, and without liability, revoke all or any portions of its acceptance of Customer's Purchase Order and/or extend any date upon which any performance thereunder is due.

8. Limited Warranty/Disclaimer of Warranties: The parties expressly agree that Customer's sole and exclusive remedy against DESIN shall be for the repair and replacement of defective Goods as provided herein. THE WARRANTIES SET FORTH BELOW ARE IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, INCLUDING THOSE OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, ALL OF WHICH ARE EXPRESSLY DISCLAIMED. NO AFFIRMATION OF DESIN, BY WORDS OR ACTION, OTHER THAN AS SET FORTH IN THIS SECTION, NOR ANY ADDITIONAL WARRANTIES, EXPRESS OR IMPLIED, MADE BY DESIN TO CONSUMER PURCHASERS OF THE GOODS SHALL CONSTITUTE A WARRANTY HEREUNDER. DESIN's warranty does not apply to any Goods that have been subjected to misuse, mishandling, misapplication, neglect (including but not limited to improper maintenance), accident, improper installation, modification (including but not limited to use of unauthorized parts or attachments), or adjustment or repair performed by anyone other than DESIN or one of DESIN's authorized agents. In no event shall DESIN be responsible for the quality or performance of any goods manufactured by other parties; such Goods shall carry only the warranty of the manufacturer.

Except as otherwise set forth below, any claim by Customer with reference to the Goods sold pursuant to the Purchase Order shall be deemed waived by the Customer unless submitted in writing to DESIN within the earlier of (i) five (5) days following the date Customer discovered, or by reasonable inspection should have discovered, any claimed breach of the foregoing warranty, or (ii) thirty (30) days following the date of shipment. Any legal dispute for breach of the foregoing warranty shall be initiated within one (1) year from the date the alleged breach was discovered or should have been discovered, whichever occurs first.

Subject to the provisions of these Standard Terms of Sale, the Limited Warranty that applies to all Goods sold is as follows:
SEE LIMITED WARRANTY.

9. Limitation of Liability: DESIN WILL IN NO EVENT BE LIABLE FOR ANY CONSEQUENTIAL, SPECIAL, INCIDENTAL, INDIRECT, PAIN AND SUFFERING, OR PUNITIVE DAMAGES TO ANY PERSON, WHETHER BASED UPON BREACH OF WARRANTY, NEGLIGENCE, STRICT LIABILITY, TORT, BREACH OF CONTRACT OR ANY OTHER THEORY, REGARDLESS OF WHETHER THE REPLACEMENT OR REFUND REMEDY SET FORTH ABOVE FAILS OF ITS ESSENTIAL PURPOSE OR FOR ANY OTHER REASON WHATSOEVER. "CONSEQUENTIAL DAMAGES" FOR PURPOSES HEREOF SHALL INCLUDE, WITHOUT LIMITATION, LOSS OF USE, INCOME OR PROFIT, OR LOSSES SUSTAINED AS THE RESULT OF INJURY (INCLUDING DEATH) TO ANY PERSON OR LOSS OR DAMAGE TO PROPERTY (INCLUDING, WITHOUT LIMITATION, PROPERTY HANDLED OR PROCESSED BY THE USE OF THE GOODS).

10. Safety Information: DESIN has supplied or made available to Customer certain information relating to safety conditions and procedures concerning the Goods, and Customer acknowledges receipt thereof and further acknowledges that it has reviewed and understands same. Customer agrees to abide by all safety procedures and warnings provided by DESIN, and further agrees to use the Goods only in conformity with the instructions provided and the intended use(s) of the Goods.

11. Trademarks: The purchase of Goods from DESIN does not entitle Customer to use, register, or otherwise identify Customer or its business with the name, trademark, service mark or other identity of DESIN or the name Obi or its likeness without express written permission from DESIN. All such marks and goodwill associated with such marks remain the sole and exclusive property of DESIN.

12. Taxes: Unless otherwise specifically provided on the face of any Purchase Order, the price for the Goods purchased does not include sales, use, excise or similar taxes, whether tariffs, duties, special assessments. The amount of any such taxes, tariffs, duties and special assessments applicable to the Goods shall be paid by Customer in the same manner and with the same effect as if originally included in the purchase price.

13. Additional Charges: If substitute or additional Goods, or repair parts, are purchased by Customer from DESIN, the provisions of these Standard Terms of Sale shall be applicable to such purchase, the same as if such substitute or additional Goods or repair parts had been originally purchased hereunder.

14. Entire Agreement: These Standard Terms of Sale, together with the terms and conditions contained on the face side of DESIN's invoice or quotation, constitute the entire agreement between DESIN and Customer with respect to the matters contained therein, and

supersedes all prior oral or written representations, proposals, correspondence, discussions, negotiations and agreements. No course of prior dealings, acceptance or acquiescence in a course of performance and no usage of the trade shall be relevant to supplement, explain or modify any terms contained herein. All representations, promises, warranties or statements by an agent or employee of DESIN that differ in any way from the terms and conditions hereof shall be given no effect or force. No waiver or alteration of terms herein, including any added or different terms identified on any invoice, quotation, or Purchase Order, shall be binding unless in writing signed by an executive officer of DESIN.

15. Priority: In the event of any inconsistency among the provisions in these Standard Terms of Sale or the Purchase Order, precedence shall be given first to these Standard Terms of Sale, and second, to the special terms and conditions contained on the face of the Purchase Order and accepted, in writing, by DESIN.

16. Governing Law: The sale of Goods pursuant to the Purchase Order shall be subject to and construed in accordance with the laws of the State of Michigan, United States of America, applicable to contracts made and performed in the State of Michigan without giving effect to any choice or conflict of law provision or rule (whether of the State of Michigan or any other jurisdiction) that would cause the application of the laws of any jurisdiction other than the State of Michigan. Exclusive jurisdiction and venue for any action arising between the parties with respect to the Goods shall lie in the Courts of Oakland County, Michigan or in federal courts situated in the Eastern District of Michigan (Southern Division) if federal court jurisdiction is invoked. The parties irrevocably consent and submit to the personal jurisdiction of said courts and agree not to challenge or assert any defense to the exclusive jurisdiction and venue of said courts including, without limitation, forum non conveniens. Customer and DESIN agree that, notwithstanding these exclusive jurisdiction and venue provisions, the parties intend to be bound by the mandatory arbitration clause articulated in these Standard Terms of Sale. The parties hereby exclude the application of the United Nations Convention on Contracts for the International Sale of Goods.

17. Severability: If any provision of these Standard Terms of Sale is deemed invalid or unenforceable by any court of competent jurisdiction or under any statute, regulation, ordinance, executive agreement or other rule of law, such provision shall be deleted or modified, at the election of Customer and DESIN, but only to the extent necessary to comply with such ruling, statute, regulation, ordinance, agreement or rule, and the remaining provisions of these Standard Terms of Sale and the Purchase Order shall remain in full force and effect.

18. Non-Waiver: No change, modification or waiver of any provision of these Standard Terms of Sale shall be valid or binding unless it is accepted in writing, by DESIN. A waiver by either party of any breach or failure to enforce any term or condition of these Standard Terms of Sale shall not in any way affect, limit or waive such party's right at any time to enforce strict compliance with that or any other term or condition of these Standard Terms of Sale.

19. Binding / Non-Assignment: This contract shall be binding upon and shall inure to the benefit of the successors, and permitted assigns of Customer and DESIN, provided, however, that Customer may not assign any portion of its rights or delegate any portion of its obligations hereunder or under any Purchase Order without in each instance obtaining DESIN's prior written consent. Any attempted assignment by Customer or assignment by Customer of this contract without DESIN's prior written consent is void.


20. DESIN's Performance Excused Upon Customer's Breach: In addition to the rights and remedies conferred upon DESIN by law or in equity, DESIN shall not be required to proceed with the fulfillment of any Purchase Order or the performance of any contract if Customer is in default in the performance of any order or contract with DESIN, and in cases of doubt as to Customer's financial responsibility, shipments under this order may be suspended or sent sight draft with bill of lading attached by DESIN. Additionally, in the case of subsequent orders or, in the case of open Purchase Orders, DESIN is not required to proceed with the fulfillment of any subsequent order or perform such open Purchase Order if Customer is in breach with regard to any prior fulfilled orders, including, without limitation, non-payment.

21. Binding Arbitration: Any controversy or claim arising out of or relating to the Goods, these Standard Terms of Sale, or any related agreement between Customer and DESIN shall be settled by binding arbitration in accordance with the following provisions:

- a. The agreement of the parties to arbitrate covers all disputes of every kind relating to or arising out of these Standard Terms of Sale and the goods, as well as any claim based on tort or any other causes of action relating to the transactions contemplated by of these Standard Terms of Sale, such as claims based on an allegation of fraud or misrepresentation and claims based on a federal or state statute. In addition, the arbitrators selected according to procedures set forth below shall determine the arbitrability of any matter brought to them, and their decision shall be final and binding on the parties.
- b. The forum for the arbitration shall be Oakland County, Michigan (USA).
- c. The governing law for the arbitration shall be the laws of the State of Michigan, as more fully provided for in the "Governing Law" section of these Standard Terms of Sale.
- d. There shall be three arbitrators, unless the parties are able to agree on a single arbitrator. In the absence of such an agreement within 10 days after the initiation of an arbitration proceeding, each party shall select one arbitrator, and those two arbitrators shall then select, within 10 days, a third arbitrator. If those two arbitrators are unable to select a third arbitrator within that 10-day period, a third arbitrator shall be appointed through the Arbitrator Selection Process facilitated through the American Arbitration Association. The decision in writing of at least two of the three arbitrators shall be final and binding on the parties.
- e. The rules of arbitration shall be the Consumer Arbitration Rules of the American Arbitration Association, as modified by any other instructions that the parties may agree on at the time, except that each party shall have the right to conduct discovery in any manner and to the extent authorized by the Federal Rules of Civil Procedure as interpreted by the federal courts. If there is any conflict between those rules and the provisions of this Section 21, the provisions of this Section shall prevail.
- f. The arbitrators' decision shall provide a reasoned basis for the resolution of each dispute and for any award. The arbitrators shall not have the power to award damages in connection with any dispute in excess of actual compensatory damages and shall not multiply actual damages or award consequential or punitive damages or award, pain and suffering, or any other damages that are excluded under the provisions of these Standard Terms of Sale.
- g. Each party shall bear its own fees and expenses (including attorney fees) regarding the arbitration and any related proceeding, and the parties shall share equally the fees and expenses of the arbitrators.
- h. The arbitrators shall have power and authority to award any remedy or judgment that could be awarded by a court of law in the state of Michigan consistent with these Standard Terms of Sale. The award rendered by arbitration shall be final and binding on the parties, and judgment on the award may be entered in any court of competent jurisdiction in the United States.



Designed, manufactured and distributed by:



7018 AC Skinner Parkway, Suite 270
Jacksonville, FL 32256
1-844-HelpObi
service@meetobi.com
www.meetobi.com